


# LUNG SPHEROID FLOW CYTOMETRY



## Purpose

The purpose of SOP 1.0 is to describe how to analyze cell subsets using flow cytometry (FC)

## Scope

SOP 1.0 is intended to cover all resources, personnel and equipment in the lung cancer organoids laboratory.

## Introduction

Flow cytometry can be used to determine different cell subsets and expression profiles. This analysis allows a high throughput characterization of cells. The unit has performed routine flow cytometry and establishment of different extra- and intracellular markers.

## 1. Cell culture media, reagents and solutions

- Phosphate-buffered saline (PBS) (Corning, New York, USA)
- Liberase™ (Sigma Aldrich)
- Flow cytometry antibodies for specific antigens (ThermoFisherScientific, Miltenyi, BD, Biolegend)
- Foxp3 / Transcription Factor Staining Buffer Set (eBioscience™, Waltham, MA USA)
- DMEM (Corning, New York, United States)

## 2.

### Equipment

- FACS Symphony (BD Biosciences)
- FACS Tubes (5ml, BD Biosciences)
- High-speed centrifuge (Hettich Rottina 35R)
- Pipette aid, serological pipettes
- FlowJo (BD Biosciences)

## 3.

### Reagent setup

- **Digestion solution:** Liberase™ (Sigma Aldrich)/Stempro (without supplements) solution 1:10
- **FACS staining buffer:** 500 ml PBS containing 1% FBS and 2mM EDTA
- **Stop solution:** DMEM containing 50 % of FBS

## 4.

### Procedure

#### 4.1

#### Human lung dissociation and generating 3D spheroids

- Single cell suspension from lung cancer spheroids: Transfer 20-40 (for larger FACS panels a greater number of spheroids are needed) spheroids in a 5ml FACS tube and incubate the tissue 15min at 37°C on a shaker with 1ml Liberase™ (Sigma Aldrich) in PBS 1:10.
- Gently pipette cells up and down to loosen the cells
- Add 2 ml of stop solution
- Centrifuge at 300 RPM for 5 minutes
- Decant the supernatant carefully
- Wash cells with 2ml PBS containing 1% FBS and 2mM EDTA
- Centrifuge at 300 RPM for 5 minutes
- Decant the supernatant carefully

- Repeat wash steps two times with 2ml of FACS staining buffer
- Decant the supernatant carefully
- Bring to a final volume to 100 µl with FACS staining buffer
- Add pre-titrated volume of selected antibodies to the suspension (use standard dilutions as recommended by reagent providers)
- Add 2µl of 100% animal (e.g mouse) serum to inhibit unspecific binding of antibodies
- Incubate for 30 min at 4°C in the dark
- Wash cells two times with 2ml of FACS staining buffer
- Immediately measure cells with the FACS Symphony
- Analyze cells with FlowJo (See [Figure 1](#))

## 4.2

### Intracellular staining

- For intracellular staining after extracellular staining wash cells with 2ml PBS without EDTA or FBS
- Add 1ml Permeabilization/Fixation Buffer (FoxP3 Staining buffer)
- Incubate for 30 min at 4°C in the dark
- Wash two times with 2ml permeabilization buffer (FoxP3 Staining buffer)
- Bring cell into suspension in 100 µl with permeabilization buffer
- Add pre-titrated volume of intracellular FACS antibodies and 2µl of 100% mouse serum; use dilutions as recommended by reagent providers.
- Incubate for 30 min at 4°C in the dark
- Wash two times with 2ml Permeabilization Buffer (FoxP3 Staining buffer)
- Immediately measure cells with the FACS Symphony
- Analyze cells with FlowJo

# 5.

## Figures

**Figure 1.** Representative image of gating strategy in FlowJo to discriminate between tumor cells and stromal compartments

